

2. Incubate the phage and the bacteria at 37°C for 15 min to allow the phage to attach to the cells.
3. Add 2–3 mL of NZY top agarose (48°C) containing IPTG and X-gal (*see Note 10*). Plate onto NZY agar plates and place them upside down in a 37°C incubator.
4. Plaques should be visible after 6–8 h. Background plaques are blue, recombinant plaques are clear and should be 10- to 100-fold above the background.
5. Count the plaques and calculate the titer. Primary libraries can be unstable. Immediate amplification of at least a portion of the library is recommended to produce a large, stable quantity of a high-titer stock of the library.

3.11. Amplification of the Primary Library

After amplification, the library is suitable for screening by a variety of techniques (2). More than one round of amplification is not recommended, since slower growing clones may be significantly underrepresented.

1. Prepare the host strains (*see Subheading 3.10.1.*).
2. Mix aliquots of the packaged library containing ~50,000 plaque-forming units (PFU) ($\leq 300 \mu\text{L}$ vol) with 600 μL of host cells in Falcon 2059 polypropylene tubes. Usually, 1×10^6 PFU are amplified (20 tubes).
3. Incubate the tubes containing the phage and host cells for 15 min at 37°C.
4. Mix 8.0 mL of melted NZY top agarose (48°C) with each aliquot of infected bacteria and spread evenly onto a freshly poured 150-mm NZY plate.
5. Incubate the plates at 37°C for 6–8 h. Do not allow the plaques to grow larger than 1–2 mm.
6. Overlay the plates with 8–10 mL of SM buffer. Store the plates at 4°C overnight with gentle rocking. The phage will diffuse into the SM buffer.
7. Recover the SM buffer containing the bacteriophage from each plate and pool it in a sterile polypropylene container. Add chloroform to a 5% final concentration and mix well.
8. Incubate for 15 min at room temperature.
9. Remove the cell debris by centrifugation for 10 min at 500g.
10. Recover the supernatant and transfer it to a sterile polypropylene container. Add chloroform to a 0.3% final concentration and store at 4°C.
11. Check the titer of the amplified library by making serial dilutions in SM buffer and plating on host cells (*see Subheading 3.10.3.*). The average titer is usually 10^9 – 10^{12} PFU/mL.
12. Frozen stocks can be made by adding dimethylsulfoxide (DMSO) to a final concentration of 7%, mixing well, and freezing at –80°C.

4. Notes

1. The mRNA sample must be highly purified for efficient cDNA synthesis. The mRNA sample may contain inhibitors that can be removed by phenol:chloroform extractions. The presence of DNA or rRNA will give an inaccurate concentration of mRNA, leading to an insufficient amount of sample used. Treat the mRNA