

Table 1
Preparation of Sucrose Solutions for 6–24% (w/v) Sucrose Gradients

	% (w/v) Sucrose solutions			
	6%	12%	18%	24%
60% (w/v) Sucrose	2.5 mL	10 mL	15 mL	20 mL
200 mM Potassium phosphate buffer, pH 6.8	12.5 mL	25 mL	25 mL	25 mL
20% (v/v) Triton X-100	125 μ L	250 μ L	250 μ L	250 μ L
H ₂ O	9.875 mL	14.75 mL	9.75 mL	4.75 mL

2. Materials

1. Borate buffer: 0.5M boric acid (H₃BO₃); bring to pH 9.0 with NaOH. Store at room temperature indefinitely. Chill the buffer to 4°C prior to virus purification.
2. 20% (v/v) Triton X-100: Add 20 mL of Triton X-100 to 80 mL of H₂O and autoclave for 5 min. When the solution cools to approx 60°C, swirl intermittently to prevent two phases from forming. Store at 4°C indefinitely.
3. Potassium phosphate buffer: Mix equal quantities of 200 mM solutions of KH₂PO₄ and K₂HPO₄ and store at 4°C. The pH of this solution should be 6.8. Dilute 200 mM potassium phosphate to obtain 50 mM and 10 mM potassium phosphate buffers and store at 4°C.
4. 60% (w/v) Sucrose: Dissolve 60 g sucrose to 100 mL with H₂O; store frozen or at 4°C.
5. 20% (w/v) Sucrose pad: Mix 30 mL of 60% (w/v) sucrose, 54 mL of 0.5M borate buffer, pH 9.0, 4.5 mL of 20% (v/v) Triton X-100 and 1.5 mL of H₂O.
6. 6–24% (w/v) Sucrose gradients: Prepare 24, 18, 12, and 6% (w/v) sucrose solutions (*see Table 1*).
7. Proteinase K (2 mg/mL) in H₂O. Store at –20°C.
8. Bentonite, prepared as described by Fraenkel-Conrat et al. (14): Stir 5 g of bentonite in 100 mL of H₂O, making sure that the bentonite is suspended well. Centrifuge at 600g for 10 min. Recover the supernatant and recentrifuge at 8000g for 20 min. Resuspend the pellet in 100 mL of 100 mM disodium ethylenediaminetetraacetate (EDTA), pH 8.0. Break the pellet up well with a glass rod and stir for 48 h at room temperature. Perform the two centrifuge steps as described previously. Resuspend the pellet obtained following the second centrifuge step in 100 mL of 10 mM sodium acetate, pH 6.0. Break the pellet up and stir overnight at room temperature. Centrifuge at 8000g for 20 min. Resuspend the pellet in 10 mL of 10 mM sodium acetate, pH 6.0. Determine the concentration of bentonite by adding 1.0 mL of bentonite to a preweighed weighing boat. Dry off the liquid in an oven at 60°C overnight and determine the weight of the bentonite. Dilute the bentonite to 10 mg/mL with 10 mM sodium acetate, pH 6.0. Store the bentonite stock at –20°C.