

The tobroviral CP is the 5' proximal gene on RNA-2, but it is not translated efficiently from the genomic RNA and is expressed from a 3' coterminal subgenomic RNA (4,5). The CP is thought to be an important determinant of nematode transmission (14). Additional ORFs have been identified on RNA-2 of several tobrovirus isolates. Some TRV RNA-2 species have undergone recombination in the 3' region with part of the same region of RNA-1, and thus carry part or all of the RNA-1-coded 29K and 16K genes (3–5,15). TRV isolate TCM and PEBV isolate SP5, RNA-2, contain an ORF for a 29.1- and a 29.6-kDa protein, respectively (5,8). Sequence comparison of RNA-2 of nematode-transmittable and nematode-nontransmittable isolates of PEBV suggest that the 29.6-kDa protein may be involved in nematode transmission (16). There is evidence that the additional 9- and 23-kDa ORFs on PEBV (SP5) RNA-2 may also have some function in vector transmission (S. A. MacFarlane, personal communication).

2. Materials

2.1. Purification of Virus

1. Celite.
2. 30 mM Sodium phosphate buffer, pH 7.5.
3. Virus extraction buffer: 50 mM sodium phosphate buffer, 0.15% thioglycollic acid (v/v). Add 1.5 mL of thioglycollic acid to 250 mL 200 mM Na₂HPO₄ and 550 mL of distilled H₂O, adjust to pH 7.5 with NaOH, and make up to 1 L with distilled H₂O. Make buffer up fresh and keep at 4°C. (**Caution:** Wear gloves when handling thioglycollic acid and use in a fume hood.)
4. Commercial blender.
5. Muslin.
6. Polyethylene glycol (PEG); mol wt = 6000.
7. Sodium chloride (NaCl).
8. Chloroform:butan-1-ol 1:1 (v/v). (**Caution:** Wear gloves and use solvents in a fume hood.)
9. Sucrose cushion: 30% sucrose in 30 mM sodium phosphate buffer, pH 7.5.
10. 10 mM Tris-HCl, pH 7.5.

2.2. Extraction of Viral RNA

1. Deoxyribonuclease I (DNase I) RNase-free.
2. 1X DNase reaction buffer: 50 mM Tris-HCl, pH 7.9, 5 mM MgCl₂.
3. EDTA.
4. Proteinase K.
5. 1X Proteinase K buffer: 100 mM Tris-HCl, pH 7.4, 1 mM EDTA, 0.1% SDS.
6. Tris-equilibrated phenol, pH 8.0 (**Caution:** Wear gloves and use in fume hood).
7. Chloroform/isoamyl alcohol 24:1 (**Caution:** Wear gloves and use in fume hood).
8. Phenol/chloroform:isoamyl alcohol 1:1.