

9. Precipitate the viral ssDNA by adding 2.5 vol of ice-cold ethanol. After at least 30 min at -20°C , pellet the precipitated nucleic acid in a microcentrifuge at 4°C . Wash the pellet with 70% ethanol, then air-dry at room temperature. Resuspend the ssDNA in sterile TE buffer.

3.2. Extraction of Total Nucleic Acids from Infected Plants

1. Grind 5–10 g of infected leaf material to a fine powder in liquid nitrogen with a mortar and pestle (*see Note 3*).
2. Suspend the frozen powdered leaf material in an equal volume of DNA extraction buffer and stir until the mixture reaches room temperature (*see Note 4*). If the homogenized plant material is too viscous, add a little more extraction buffer.
3. Centrifuge the mixture for 10 min at 4°C in 30-mL polypropylene tubes at 12,000g (in a Sorvall SS34 fixed angle rotor) to pellet the plant debris.
4. Transfer the supernatant to 30-mL polypropylene tubes containing 10 mL of Tris-buffered phenol and mix well. Centrifuge at 4°C for 10 min at 12,000g.
5. Transfer the aqueous phase into clean polypropylene tubes (*see Note 5*). Add an equal volume of chloroform, mix, and centrifuge again at 4°C for 10 min at 12,000g. Repeat this step.
6. Transfer the aqueous phase to clean polypropylene tubes. Add 0.7 vol of isopropanol and mix gently. A stringy precipitate of nucleic acids should appear in the tube.
7. Pellet the nucleic acids at 17,000g for 20 min, discard the supernatant, and stand the tubes upside down to drain off excess isopropanol.
8. Wash the pellet by adding about 10 mL of 70% ethanol, dislodge the pellet, and mix and centrifuge again at 17,000g for 5 min at 4°C (*see Note 6*).
9. Invert the centrifuge tube on absorbent paper and allow the pellet to air-dry for about 10 min at room temperature.
10. If only a crude nucleic acid preparation is required, resuspend the pellet in 0.5–1 mL of TE buffer (*see Note 7*).

3.3. Isolation of RF-DNA by Alkaline Denaturation and Anion-Exchange Chromatography

This procedure uses a plasmid isolation protocol based on alkaline denaturation of chromosomal DNA and purification of cccDNA by anion-exchange chromatography. The protocol is basically as described by the manufacturers of the kit that we use (the Qiagen-tip 20 kit, Qiagen GmbH and Qiagen), but can easily be adapted for the use of similar plasmid isolation kits.

1. Resuspend the total nucleic acid pellet from **Subheading 3.1., step 9** in 0.3–0.5 mL of solution P1 in a microcentrifuge tube (*see Note 8*).
2. Add the same volume of solution P2. Mix gently to avoid shearing chromosomal DNA and incubate at room temperature for 5 min.
3. Add the same volume (0.3–0.5 mL) of ice-cold solution P3. Mix well and place on ice for 10 min.