

2. 10X DIG-DNA labeling mixture: 1 mM dATP, 1 mM dCTP, 0.65 mM dTTP, 0.35 mM Dig-dUTP, pH 7.5 (20°C) (Boehringer Mannheim). Store at -20°C.
3. Klenow enzyme (2 U/ μ L) (Boehringer Mannheim). Store at -20°C.
4. 4M LiCl. Store at 4°C.

2.3. Dot-Blot Hybridization

1. Mortars.
2. 50 mM Sodium citrate, pH 8.5. Autoclaved. Store at room temperature.
3. Nylon membrane, positively charged (Boehringer Mannheim).
4. 20X SSC: 3M NaCl, 0.3M sodium citrate, pH 7.0. Autoclaved. Store at room temperature.
5. Blocking stock solution: 10% (w/v) blocking reagent (Boehringer Mannheim) in buffer 1 (*see Subheading 2.4., item 1*); autoclaved and stored at 4°C.
6. Prehybridization solution: 50% (v/v) deionized formamide, 5X SSC, 2% (w/v) blocking reagent (added from 10% blocking stock solution), 0.1% (w/v) *N*-lauroylsarcosine, 0.02% (w/v) sodium dodecyl sulfate (SDS). Store at -20°C.
7. Washing solution I: 2X SSC, 0.1% (w/v) SDS.
8. Washing solution II: 0.1% SSC, 0.1% (w/v) SDS.

2.4. Immunological Detection

1. Buffer 1: 0.1M maleic acid, 0.15M NaCl, pH 7.5 (20°C), adjusted with concentrated NaOH, autoclaved.
2. Washing buffer: 0.3% Tween-20 (Sigma, St. Louis, MO) in buffer 1.
3. Buffer 2: 1% blocking reagent (added from 10% sterile stock solution) in buffer 1.
4. Buffer 3: 0.1M Tris-HCl, pH 9.5, 0.1M NaCl, 50 mM MgCl₂. Discard the solution if a precipitate appears after a long storage period.
5. Antibody-conjugated solution freshly prepared: Dilute 1:10,000 the sheep antidigoxigenin Fab-fragments conjugated to alkaline phosphatase (750 U/mL) (Boehringer Mannheim), in buffer 2. Diluted solution is stable only for 12 h at 4°C. Undiluted preparation are stable at 4°C. Do not freeze.
6. 75 mg/mL Nitroblue tetrazolium (NBT), in 70% (v/v) dimethylformamide, and 50 mg/mL 5-bromo-4-chloro-3-indolyl phosphate (BCIP), in dimethylformamide (Boehringer Mannheim). Store at -20°C.
7. CSPD® substrate (10 mg/mL) (Tropix, Bedford, MA). Store at 4°C and protected from light.
8. X-Ray films.

3. Method

3.1. Synthesis of the Probe

The two methods described here for preparing the nonisotopic label require that part of the viral sequence is cloned into the polylinker site of a plasmid vector that contains a promoter for SP6, T7, or T3 RNA polymerase. For RNA probes, the plasmid vector must be digested with an appropriate restriction