

3. Resuspend the pellet by gentle shaking in 200 μL of 0.03M sodium phosphate, pH 8.0, per 100 g leaf tissue, and then centrifuge for 10 min at 12,000g in a cooled Eppendorf centrifuge. Using a Pasteur pipet, carefully transfer the virus-containing supernatant to a new 1.5-mL centrifuge tube. The preparation at this step, although containing still a little host-plant contamination, is good for viral RNA isolation and subsequent cDNA cloning. The quality can be controlled in an electron microscope and spectrophotometrically. The $\text{OD}_{260/280}$ of the preparation should range from 1.45 to 1.60. The concentration can be estimated from the extinction coefficient $E_{260\text{ nm}} \approx 5.1\text{--}5.3$, which is equivalent to 1 mg/mL at 1 cm length of the light path.

3.5. Sucrose Gradient Centrifugation

If the preparation is to be used for antiserum production, a further purification step by sucrose density gradient centrifugation is recommended. Load at maximum 1 mL, with up to 3 mg of the virus suspension per each preformed sucrose gradient, ranging from 10 to 40% sucrose in 0.03M sodium phosphate, pH 8.0, and centrifuge for 4 h at 113,000g in a Beckman SW-28 rotor. Carefully collect the virus-containing zones from the gradients with an ISCO density gradient fractionator, according to the absorption profile, monitored at 254 nm. Dilute the virus-containing fractions four times in 0.03M sodium phosphate, pH 8.0, and sediment the virus by centrifugation at 123,000g for 2.5 h in a Beckman Ti-60 rotor. Resuspend the final pellet in 0.2 mL of the same buffer.

3.6. Storage of Purified Virus

For short periods, i.e., a few days, purified virus is stored at 4°C at the highest concentrations possible. For long-term storage, the purified virus should be aliquoted in useful amounts and stored frozen at -80°C or in liquid nitrogen.

3.7. Extraction of Viral RNA (see Note 11)

1. Mix the virus preparation with 1 vol of twofold proteinase K buffer and add proteinase K to an end concentration of 400 $\mu\text{g/mL}$. Incubate the mixture at 37°C for 30 min to digest the viral coat protein.
2. Add 1 vol of TE phenol to the digestion mixture and strongly vortex for 30 s; break the emulsion by a centrifugation for 4 min at room temperature and 8000g in a table centrifuge equipped with a swingout rotor. Using a pipet, carefully remove the upper aqueous phase, which contains the RNA.
3. Extract the aqueous phase again, first with phenol:chloroform mix, and then with the chloroform mix, as described above. Each time, carefully remove the aqueous phase.
4. Add 0.1 vol 3M sodium acetate, pH 5.2, and 2.5 vol of ice-cold 100% ethanol to aqueous phase and incubate for at least 30 min at -70°C. This will precipitate the RNA from the solution.