

**Note 2).** The solution is gently stirred until the PEG is completely dissolved, then held without stirring for 1 h at 4°C. After centrifugation for 30 min at 12,000g, the pellets are resuspended in a maximum of 3 mL of buffer C (*see Note 3*).

### 3.2.4. Separation on Sucrose Gradient

Sucrose gradients 10–40% are prepared from two stock sucrose solutions (*see Subheading 2.1., items 4 and 5*). Viral pellets, resuspended in buffer C, are loaded directly onto six gradients (500  $\mu$ L per gradient) and ultracentrifuged for 2.5 h at 300,000g, in a Beckman SW 41 rotor. Gradient tubes are scanned at 254 nm with a ISCO ultraviolet absorption monitor and the UV absorbing regions corresponding to virus fractions are collected in ultracentrifuge tubes.

### 3.2.5. Concentration

To concentrate the virus, the virus-containing fractions are ultracentrifuged 16 h at 95,000g in a Beckman 60 Ti rotor. The final virus pellet is then resuspended in 200  $\mu$ L of buffer C. The virus concentration can be estimated by measuring the absorption at 254 nm. Concentration is then calculated based on the formula: 2.4 U of OD<sub>254nm</sub> correspond to a concentration of 1 mg/mL of virus.

## 3.3. Viral RNA Isolation

### 3.3.1. RNA Extraction

The virus suspension is incubated in the presence of proteinase K (200 ng/ $\mu$ L) and 0.5% SDS for 15 min at 50°C. This suspension is then extracted with an equal volume of Tris-HCl saturated-phenol:chloroform:isoamyl alcohol (25:24:1, [v/v/v]), and then with an equal volume of chloroform:isoamyl alcohol (24:1, [v/v]). The aqueous phase is ethanol-precipitated in the presence of 0.1 vol of 3M sodium acetate, pH 5.3, and 2.5 vol of 96% ethanol at –20°C, for at least 1 h. Viral RNA is recovered by centrifugation for 20 min at 4°C at 15,000g. The RNA pellet is finally resuspended in 10–20  $\mu$ L of TE buffer or water.

### 3.3.2. Analysis of Viral Nucleic Acid

RNA is denatured with formaldehyde and electrophoresed in a 0.8% agarose formaldehyde gel, as described by Miller (21), with few modifications (*see Note 4*). The RNA volume should be reduced to 1  $\mu$ L, to which 2.2  $\mu$ L of buffer D and 5.8  $\mu$ L of buffer E (*see Subheading 2.3., item 3*) are added. After denaturation by heating at 70°C for 10 min, 1  $\mu$ L of gel-loading buffer is added and the sample is loaded on the gel (*see Subheading 2.3., items 5–7*). After electrophoresis, the gel is stained for 15 min and destained in water (*see Subheading 2.3., item 8*).