

3. Allow all components, except the burst disks, to soak for at least 1 h to sterilize. Sterilize the burst disks by soaking for 5 min. After sterilization, transfer all components to a layer of autoclaved wipes to dry.

3.2.1. Preparation of Microcarriers

The following method is modified slightly from that of Heiser (5) and was originally developed by Sanford et al. (4). It allows the preparation of 30 mg of particles, which would be adequate for 50 individual bombardments. We usually prepare 30 mg, since the particles can be stored for only about 3 wk at 4°C, or 2 wk at room temperature. For ease of handling, when 60 mg are needed, we prepare two tubes with 30 mg each.

1. In a 1.5-mL microcentrifuge tube, vigorously vortex 30 mg of particles in 0.5 mL freshly prepared 70% EtOH for 5 min.
2. Incubate particles at room temperature for 15 min.
3. Microcentrifuge (15,000g) for 5 s and decant.
4. Wash particles three times as follows:
 - a. Add 0.5 mL sterile water.
 - b. Vortex for 1 min.
 - c. Allow particles to settle for 1 min.
 - d. Microcentrifuge for 5 s and decant.
5. Add sterile 50% glycerol to bring the particle concentration to 60 mg/mL.
6. While vigorously vortexing, distribute 40- μ L aliquots of particles into 1.5-mL microcentrifuge tubes.

3.2.2. Coating DNA on Microcarriers

1. Add 1 and 3 μ L (1 mg/mL) CsCl₂-purified plasmid encoding the bialaphos-resistance gene and the gene of interest, respectively, to a 40 μ L aliquot of particles (see Note 2). We usually bombard with four different constructs in a given experiment, which results in coating four aliquots of particles at a time.
2. Place tubes on a table mixer (e.g., Eppendorf Model 5432) and mix for 5 min.
3. While vortexing vigorously, add 40 μ L CaCl₂ (2.5M) and then 16 μ L spermidine (0.1M).
4. Continue vigorous vortexing for 3 min.
5. Allow particles to settle (for at least 1 min) while subsequent aliquots are being processed.
6. Microcentrifuge for 4–5 s.
7. Draw off supernatant with a pipet and discard.
8. Without disturbing pellet, add 112 μ L freshly prepared 70% EtOH.
9. Draw off supernatant with a pipet and discard.
10. Again without disturbing pellet, add 112 μ L 100% EtOH.
11. Draw off supernatant with a pipet and discard.
12. Suspend particles in 39 μ L of 100% EtOH (a few microliters will be lost to evaporation and 36 μ L are needed) by mixing with the pipet tip and then vigorously pipeting up and down.