

2. Centrifuge at 10,000g for 20 min at 4°C.
3. Filter aqueous phase through large kimwipes supported by cheesecloth. Make sure that no chloroform or carbon tetrachloride remains in this initial extract (*see Note 3*).
4. Add PEG (mol wt 8000) to 4% (w/v) and NaCl to 1.75% (w/v). Stir on ice for 1 h (*see Note 4*).
5. Centrifuge at 10,000g for 15 min at 4°C. Carefully discard supernatant and retain the pellets.
6. Resuspend pellets in ca. 0.1 vol of initial extract, using cold 0.5M borate, pH 8.0 (*see Note 5*).
7. Just before low-speed centrifugation, add Triton X-100 to 0.25% (v/v) and stir until thoroughly dissolved. Low-speed centrifuge at 8000g for 15 min at 4°C. Retain the supernatant.
8. Ultracentrifuge the retained supernatant at 70,000g for 1.5 h at 4°C (*see Note 6*).
9. Resuspend pellets in buffer C, using about one-thirtieth vol of the initial extract (*see Note 5*).
10. Repeat **steps 7–9**, but resuspend pellets in 5 mL buffer C. Centrifuge briefly at low speed (8000–10,000g) for 1–2 min.
11. Layer supernatant on 10–40% sucrose density gradients made in buffer C, and centrifuge 2 h at 95,000g in a SW 28 swing bucket rotor (*see Notes 7–9*).
12. If desired, dilute virus-containing fractions with at least 3 vol of buffer C and mix thoroughly. Ultracentrifuge as in **step 8**, above, to concentrate virus-containing fractions.
13. Resuspend pellet in 1 or 2 mL of buffer C (*see Note 10*).

### 3.2. Purification of Maize Dwarf Mosaic Virus

This method, originally developed by W. Langenberg, has been used by a number of laboratories for purification of potyviruses infecting Gramineae, particularly maize dwarf mosaic (MDMV) and the related sugarcane mosaic virus (SCMV), sorghum mosaic virus (SrMV), and johnsongrass mosaic virus (JGMV). Langenberg reports that addition of 2M guanidine HCl, at **steps 1–6**, will significantly increase virus yield (personal communication).

1. Collect 500–1500 g infected tissue and grind in a large, prechilled blender with a minimum of 1 L buffer B. Strain the contents through several layers of cheesecloth, squeezing out as much sap as possible (*see Notes 11–13*).
2. Add carbon tetrachloride to 5% (v/v) and mix in a blender for about 5 s.
3. Centrifuge for 10 min at 10,000–15,000g at 4°C and reserve the supernatant (*see Note 14*).
4. While stirring, add Triton X-100 to 0.25% final concentration. Then, add solid PEG to 6% final concentration. Stir until dissolved (about 0.5 h) (*see Notes 15 and 16*).
5. Centrifuge solution for 20 min at 10,000–15,000g at 4°C and discard supernatant. Resuspend pellet in ca. 100 mL buffer C and clarify by centrifuging for 10 min at 10,000g at 4°C, reserving the supernatant.