

Nelson and Tremaine (8) developed an alternative protocol in which the tissue was homogenized in an acidic sodium acetate buffer to achieve clarification of the extract in the absence of organic solvents. A later improvement of this method included a polyethylene glycol precipitation step (9) and has been widely used for carmovirus purification (*see Note 6*). Here, we describe a different modification of the acetate buffer method that renders high yields of high-quality virus in a significantly shorter time. The method described here is based on a high-speed centrifugation to pellet the virions through a sucrose cushion. Our experience with several carmoviruses shows that the virions obtained in this way are pure enough for different virological and molecular analyses (*see Notes 4 and 5*). Viral RNA extraction from purified virions is carried out essentially as described previously (10).

2. Materials

All the solutions should be prepared with MilliQ water.

2.1. Virus Purification

1. Liquid nitrogen.
2. 0.2M Sodium acetate, pH 5.0. Autoclave and store at room temperature.
3. 10 mM Tris-HCl, pH 7.3. Autoclave and store at room temperature.
4. 20% (w/v) Sucrose in 10 mM Tris-HCl, pH 7.3. Freshly prepared in each experiment. Alternatively, the solution can be autoclaved at 121°C for no more than 15 min to avoid caramelization of sucrose.

2.2. Viral RNA Extraction

1. 4 mg/mL Proteinase K dissolved in sterile water. Store aliquoted at -20°C.
2. 10% (w/v) Sodium dodecyl sulfate (SDS) in 50 mM Tris-HCl, pH 8.0. Autoclave and store at room temperature.
3. 1M Tris-HCl, pH 7.6. Autoclave and store at room temperature.
4. 5M NaCl. Autoclave and store at room temperature.
5. Phenol:chloroform:isoamyl alcohol (25:24:1). Prepared and stored as described by Sambrook et al. (11). **Caution:** Phenol is highly toxic and gloves must be worn when handling it.
6. Diethyl ether or chloroform. **Caution:** Diethyl ether is highly volatile and should be stored and used in a fume hood.
7. 3M Sodium acetate, pH 5.5. Prepared and stored as described (11).
8. Absolute ethanol and 70% (v/v) ethanol in water. Store at -20°C.

3. Methods

When possible, keep the extract on ice or at 4°C. Also, once the virus purification (and RNA extraction) has been started, it should be carried out as quickly as possible. **Caution:** Caution must be taken with RNase contamination; gloves