

10. Macrocarrier loader (Bio-Rad).
11. Box of laboratory wipes (Kimwipes) wrapped in aluminum foil and autoclaved.
12. Gold microcarriers (1.0 μm) (Bio-Rad).
13. Plasmid DNA encoding a chimeric phosphinothricin acetyl transferase gene (**16**) (e.g., we have used either the CaMV 35S [Huntley and Hall, in preparation] or maize ubiquitin 1 [Christensen and Quail, in preparation] promoter, coupled with the nopaline synthase polyadenylation signal, to drive expression in rice).
14. Plasmid DNA encoding the gene of interest to be introduced.
15. 0.1M Spermidine. Dissolve the free base in water, filter-sterilize, and store at -20°C . Discard after 1–2 mo.
16. CaCl_2 (2.5M), filter-sterilize, and store at -20°C , stable for months.
17. 50% Glycerol, filter-sterilize; make up as needed.

2.3. Selection and Regeneration

1. LS 2.5 tissue-culture medium supplemented with 4 mg/L bialaphos.
2. MSD4 tissue-culture medium:
 - a. To 985 mL ultrapure H_2O , add the following with stirring: 1 package Murashige and Skoog Salt Mixture ([15], Gibco-BRL, cat. no. 11117-066), 30 g sucrose, 100 mg myo-inositol, 1 mL of MS vitamins + glycine (1000 X stock, *see below*), 0.5 mL BAP (6-benzylaminopurine, Sigma cat. no. B 9395; 1 mg/mL stock in 0.1M HCl), 50 μL NAA (naphthalene acetic acid, Sigma cat. no. N 0640; 1 mg/mL stock in EtOH). Adjust pH to 5.80 with 0.5M KOH. Add 3.85 g agarose (Sigma Type I, cat. no. A6013) per 1-L media bottle, and add medium. Cover bottle lid and neck with foil, leave lid loose, and autoclave (121°C , 25 min). After autoclaving, swirl medium to evenly distribute agarose and cool to 50°C . Transfer bottles to a laminar flow hood, spray, and wipe with 70% EtOH. Pour approx 50 mL per $100 \times 25\text{-mm}$ plate. Leave plates overnight in laminar flow hood (with fan turned off) to dry.
 - b. MS vitamins + glycine (1000X stock; store at -20°C , stable for months): To 50 mL ultrapure H_2O , add the following, with stirring: 25 mg nicotinic acid, 25 mg pyridoxine HCl, 5 mg thiamine-HCl, 100 mg glycine.
3. MS0 medium: Prepare as for MSD4, except omit BAP and NAA. After adding the agarose and medium to the media bottle, fully dissolve the agarose by heating (Phytigel, Sigma cat. no. P 8169 may be substituted for agarose). Aliquot 60–65 mL per Magenta box (Sigma cat. no. V8505) and put the box lids on. Autoclave (121°C , 25 min) and transfer Magenta boxes directly to a laminar flow hood to cool. Stacking boxes helps eliminate condensation on the lids.

2.4. Plant Growth

1. Shredded pasteurized peat moss.
2. Coarse vermiculite.
3. 1-gal Plastic pots.
4. $12 \times 2\text{-in.}$ Round trays.
5. $1 \times 2\text{-mm}$ Mesh nylon screen.