

sequenced furoviruses: first, the 5'-proximal position of the CP cistron on a genomic RNA; second, the presence of an open reading frame (ORF) adjacent to the CP cistron, which is expressed by translational readthrough for BNYVV (13,14), SBWMV (6,15), and probably PMTV (10), and by ribosome scanning for PCV (7,16). In the case of BNYVV, this protein has been implicated in fungus transmission (17). Finally, all the characterized furoviruses have an ORF for a small cysteine-rich protein, which has been shown in the case of BNYVV to be involved in regulation of CP expression (18).

At present, a dozen viruses have been classified within the furovirus group (19), and, for each of them, a distinct purification procedure (or several such procedures) has been described. Because of this diversity, this chapter will only detail the technique described by Putz and Kuszala (*see ref. 20*), which is currently used in slightly modified form in this laboratory for large-scale purifications of BNYVV. However, significant differences with the purification methods for other well-studied furoviruses are mentioned in **Subheading 4**.

Generally speaking, furoviruses do not multiply well in their hosts (19). Another obstacle to their purification is a tendency of the virions to self-aggregate or to adhere to large cellular components. Aggregation is reported to be particularly severe in the case of PMTV and *Nicotiana velutina* mosaic virus (NVMV) (21–23). Consequently, yields of purified virus rarely exceed 30 mg/kg of leaves. However, the method described in this chapter has given yields of BNYVV of as much as 80 mg/kg of leaves of the local lesion host *Chenopodium quinoa*.

## 2. Materials

1. Plant host: *C. quinoa*, aged 7–8 wk.
2. Inoculum: Freshly ground leaves of infected *C. quinoa*. Leaves that have been frozen for some weeks or lyophilized (preferable) can also be used.
3. Inoculum buffer (4X): 0.2M potassium phosphate, pH 7.0.
4. Abrasive powder: Celite.
5. Sterile cotton, mortar, and plastic gloves.
6. Variable speed food blender of approx 1-L vol.
7. Extraction buffer: Sterile 0.1M sodium tetraborate, adjusted to pH 9.0 with boric acid.
8. Clarification agent: carbon tetrachloride (CCl<sub>4</sub>).
9. Miracloth or equivalent filter tissue.
10. Precipitation agent: NaCl (8 g/L of clarified sap), polyethylene glycol (PEG) 6000 (20 g/L of clarified sap).
11. 20% Sucrose cushion (w/v) prepared in 10 mM extraction buffer.
12. Low-speed centrifuge equipped with fixed-angle rotors and with tubes of capacity approx 0.5 L, 0.25 L, and 30 mL. Ultracentrifuge with fixed-angle rotor and accepting ultracentrifuge tubes of approx 25 mL.