

4. Tris-buffered phenol, prepared as described in **Subheading 2.1.2.** (**Caution:** Handle solutions containing phenol with caution, and in a fume hood: Phenol is toxic and highly corrosive.)
5. Chloroform (**Caution:** to be handled in a fume hood, since chloroform is carcinogenic).
6. Isopropanol (propan-2-ol).
7. Absolute ethanol and 70% ethanol.
8. TE buffer: 10 mM Tris-HCl, 1mM EDTA, pH 8.0, sterilized by autoclaving.

2.3. Isolation of RF-DNA by Alkaline Denaturation and Anion-Exchange Chromatography

1. Commercially available plasmid isolation kit, for example the Qiagen-tip 20 kit (Qiagen GmbH and Qiagen, Hilden, Germany), with Qiagen-tip 20 columns and solutions:
 - P1: 50 mM Tris-HCl, 10mM EDTA, pH 8.0, 100 µg Rnase A.
 - P2: 200 mM NaOH, 1% SDS.
 - P3: 3.0M KOAc, pH 5.5.
 - QBT: 750 mM NaCl, 50 mM MOPS, 15% ethanol, pH 7.0, 0.15% Triton X-100.
 - QC: 1.0M NaCl, 50 mM MOPS, 15% ethanol, pH 7.0.
 - QF: 1.25M NaCl, 50 mM Tris-HCl, 15% ethanol, pH 8.5.
2. Sterile 1.5-mL microcentrifuge tubes.
3. Sterile TE buffer: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0.

3. Methods

3.1. Isolation of Maize Streak Virus Virions and ssDNA

1. Homogenize freshly harvested leaf material (*see Note 2*) from infected plants with an equal weight/volume of 0.1M acetate buffer, pH 4.8, at room temperature.
2. Squeeze the homogenate through a layer of cheesecloth.
3. Immediately adjust the pH of the homogenate back to pH 4.8 with 10% glacial acetic acid.
4. Remove the precipitated plant components by low speed centrifugation (12,000g for 10 min in a Sorval GSA rotor).
5. Pellet the virions by ultracentrifugation (130,000g for 150 min in a Beckman Type 35 rotor).
6. Resuspend the high-speed pellet in 0.05M sodium phosphate buffer, pH 7.5, and repeat **steps 3** and **4**.
7. Resuspend the virus pellet in 0.05M sodium phosphate buffer, pH 7.5. Alternatively, if the virus is to be used for extraction of ssDNA, resuspend in DNA extraction buffer.
8. Isolate ssDNA from virions by extraction with an equal volume of phenol:chloroform (1:1). Remove the aqueous phase containing ssDNA to a fresh tube and re-extract with phenol:chloroform, if necessary. Remove residual phenol by extraction with chloroform.