

10. It is essential that bubbles are not present to assure uniform flow of buffer, and thereby uniform transfer of the DNA, to the membrane.
11. Parafilm or other hydrophobic material can be butted up to the gel on all sides to prevent overhanging paper from absorbing buffer without it going through the gel.
12. Use 5×10^6 cpm denatured probe/mL hybridization solution. To denature, dilute the volume of probe needed to at least 50 μ L with H₂O, heat denature (95–100°C, 5 min), and fast chill for 2–3 min on wet ice prior to adding to the hybridization solution.
13. If the membrane becomes dry, it is virtually impossible to remove the hybridized probe. If kept moist, ~90–100% of the probe can be removed and the blot reused. (We have reused nitrocellulose blots 6–7 times with only limited loss in sensitivity.) Strip the probe by incubating it in 200–300 mL of 30 mM NaOH, 1 mM EDTA for 15 min at room temperature, with gentle agitation. Decant and neutralize with a similar volume of Southern neutralization solution. Use a radiation monitor before and after stripping to evaluate the effectiveness of probe removal. If some radioactivity remains, expose the blot for a specific time to determine the background level. To re-use the blot, the prehybridization step can be omitted. Add the hybridization solution (without probe) and incubate at 65°C for 1–2 h before adding the probe.

Acknowledgments

Support from the Rockefeller Foundation Rice Biotechnology Program, the Texas A&M Office of Sponsored Projects, TAES Research Enhancement Program, and the Texas Advanced Technology Program (Projects 010366-108 and 010366-149) is gratefully acknowledged. We thank Kimberly Snowden for critical reading of this manuscript.

References

1. Southern, E. (1975) Detection of specific sequences among DNA fragments separated by gel electrophoresis. *J. Mol. Biol.* **98**, 503–511.
2. Taylor, B. H., Manhart, J. R., and Amasino, R. M. (1993) Isolation and characterization of plant DNAs, in *Methods in Plant Molecular Biology and Biotechnology* (Glick, B. R. and Thompson, J. E., eds.), CRC, Boca Raton, pp. 360.
3. Agudo, L. d. C., Gavidia, I., Perez-Bermudez, P., and Sugura, J. (1995) PEG precipitation, a required step for PCR amplification of DNA from wild plants of *Digitalis obscura* L. *Biotechniques* **18**, 766–768.
4. Xu, Y., Buchholz, W. G., DeRose, R. T., and Hall, T. C. (1995) Characterization of a rice gene family encoding root-specific proteins. *Plant Mol. Biol.* **27**, 237–248.
5. Matzke, M. A. and Matzke, A. J. M. (1995) How and why do plants inactivate homologous (trans)genes? *Plant Physiol.* **107**, 679–685.
6. Finnegan, J. and McElroy, D. (1994) Transgene inactivation: plants fight back! *Biotechnology* **12**, 883–888.